

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CD284 Monoclonal Antibody

单克隆抗体

产品货号: MA02574

第三版

免疫原:Fusion Protein of CD284

亚型:IgG2b

克隆:3G9A4

描述:TLR4,also named as CD284,belongs to the Toll-like receptor family. TLR4 cooperates with LY96 and CD14 to mediate the innate immune response to bacterial lipopolysaccharide (LPS). TLR4 acts via MYD88,TIRAP and TRAF6,leading to NF-kappa-B activation,cytokine secretion and the inflammatory response. Three alternatively spliced transcript variants that encode different protein isoforms have been described.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF