

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human CD137L Monoclonal Antibody

单克隆抗体

产品货号: MA02575

第三版

免疫原:Fusion Protein of TNFSF9

亚型:IgG1

克隆:1B4B1

描述:TNFSF9,also named as 4-1BBL,belongs to the tumor necrosis factor family. It is a cytokine that binds to TNFRSF9. TNFSF9 induces the proliferation of activated peripheral blood T-cells. It may have a role in activation-induced cell death (AICD) and may play a role in cognate interactions between T-cells and B-cells/macrophages.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF