



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat VDAC1 Monoclonal Antibody

单克隆抗体

产品货号: MA02581

第三版

免疫原:Fusion Protein of VDAC1/Porin

亚型:IgG3

克隆:1E2C7

描述:VDAC1,also named as VDAC,Porin 31HM,Porin 31HL and Plasmalemmal porin,belongs to the eukaryotic mitochondrial porin family. It adopts an open conformation at low or zero membrane potential and a closed conformation at potentials above 30-40 mV,to form a channel through the mitochondrial outer membrane and also the plasma membrane. Unlike other membrane transport proteins,porins are large enough to allow passive diffusion.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF