

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ZC3HAV1 Monoclonal Antibody

单克隆抗体

产品货号: MA02584

第三版

免疫原:Fusion protein of ZC3HAV1

亚型:IgG1

克隆:1G10B9

描述:the zinc-finger antiviral protein (ZAP) was originally identified from a rat cDNA library due to its conferring resistance to retroviruses. ZC3HAV1 is a CCCH-type zinc finger protein,it may primarily function to inhibit viral gene expression and induce an innate immunity to viral infection. Alternative splicing occurs at this locus and two variants,each encoding distinct isoforms.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF