

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CST3 Monoclonal Antibody

单克隆抗体

产品货号: MA02626

第三版

免疫原: Recombinant Protein of Cystatin C of CST3

亚型: IgG

克隆: 3B6

描述: As an inhibitor of cysteine proteinases, this protein is thought to serve an important physiological role as a local regulator of this enzyme activity. Known to inhibit cathepsin B, H, and L.

配方:

浓度: 1mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF