

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Ubiquitin Monoclonal Antibody

单克隆抗体

产品货号: MA02628

第三版

免疫原:Synthetic Peptide of Ubiquitin

亚型:IgG

克隆:6C4

**描述:**This gene encodes a member of a subfamily of the adenosine deaminase protein family. The encoded protein is one of two adenosine deaminases found in humans, which regulate levels of the signaling molecule, adenosine. The encoded protein is secreted from monocytes undergoing differentiation and may regulate cell proliferation and differentiation. This gene may be responsible for some of the phenotypic features associated with cat eye syndrome. Alternative splicing results in multiple transcript variants. ADA2 (Adenosine Deaminase 2) is a Protein Coding gene. Diseases associated with ADA2 include Sneddon Syndrome and Polyarteritis Nodosa, Childhood-Onset. Among its related pathways are Innate Immune System and Metabolism of proteins.

**配方:**

浓度:1mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF