

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat NBR1 Monoclonal Antibody

单克隆抗体

产品货号: MA02644

---

第三版

免疫原: Recombinant Protein of NBR1

亚型: IgG

克隆: 4C6

**描述:** The protein encoded by this gene was originally identified as an ovarian tumor antigen monitored in ovarian cancer. The encoded protein contains a B-box/coiled-coil motif, which is present in many genes with transformation potential. It functions as a specific autophagy receptor for the selective autophagic degradation of peroxisomes by forming intracellular inclusions with ubiquitylated autophagic substrates. This gene is located on a region of chromosome 17q21.1 that is in close proximity to the BRCA1 tumor suppressor gene. Alternative splicing of this gene results in multiple transcript variants. NBR1 (NBR1, Autophagy Cargo Receptor) is a Protein Coding gene. Diseases associated with NBR1 include Inclusion Body Myositis and Synucleinopathy. Among its related pathways are Mitophagy - animal. GO annotations related to this gene include ubiquitin binding and mitogen-activated protein kinase binding.

**配方:**

浓度: 1mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF