
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

COX-2 (C580) Anti-Rat Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01114

第三版

描述:Cyclooxygenase (COX) is a membrane bound enzyme responsible for the oxidation of arachidonic acid to Prostaglandin G2 (PGG2) and the subsequent reduction of PGG2 to PHG2.

抗原:Synthetic peptide for a part of Rat COX-2

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Rat COX-2 specific. Non-cross reacts with Human COX-2 and Rat COX-1.