
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Endothelin-A Receptor (A405) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01143

第三版

描述:Two kinds of receptors for Endothelin have been found, namely Endothelin-A (ET-A) receptor and Endothelin-B (ET-B) receptor. ET-A receptor is ET-1 specific and expressed in lung, aorta, atrium, liver, muscle, cerebral cortex, cerebellum, kidney, adrenal, duodenum, colon and placenta at the mRNA or cDNA level, while ET-B receptor lacks isoepitope specificity and is found in liver, lung, placenta, kidney, heart and brain.

抗原:Synthetic peptide of C-terminal of Human ET-A receptor

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Cross-reacts with Bovine and Rat.