

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

EPRS (Glutamyl-prolyl tRNA synthetase)

多克隆抗体

产品货号: PA01146

第三版

描述:EPRS is a component of ARS multicomplex. EPRS is released and phosphorylated by IFN stimulation and forms a complex with GAIT.

抗原:V573-K889 native human recombinant EPRS.

配方:Anti-glutamyl-prolyl-tRNA synthetase rabbit serum containing 0.05% sodium azide for preservative.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:This antibody reacts specifically with 155kDa glutamyl prolyl tRNA synthetase protein.