

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ERA (E. coli Ras-like protein) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01147

第三版

**描述:**ERA (Escherichia coli Ras-like protein) is an E. coli GTP binding protein that is essential for proliferation. A DNA database search suggests that homologous sequences with ERA exist in various organisms including human, mouse, Drosophila, Caenorhabditis elegans and Antirrhinum majus.

**抗原:**Synthetic peptide of middle part of human ERA (E. coli Ras-like protein) ((C)AVKDPNTQSVGNPQ)

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**Confirmed by western blotting using transfectant cells