

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Fatty Acid Synthase Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01148

第三版

描述:Animal fatty acid synthase (FAS) is a homodimer protein, which synthesizes long-chain fatty acids and is rich in liver, brain, breast and lung. However, the precise cellular localization of FAS in human tissue has not been elucidated.

抗原:Synthetic peptide for a part of Human Fatty Acid Synthase.

配方:Lyophilized product from 1 % BSA in PBS containing 0.05 % Na₃N

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Human FAS specific. Cross reactivity with rat and mice are not tested