
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Fezf1 Anti-Mouse Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01150

第三版

描述: Mouse Fezf1 (also called Fez, Zfp312-like or 3110069A13Rik) is a transcription factor, C2H2-type zinc finger protein consists of 475 amino acids, and it is expressed in the olfactory bulb (OB) and hypothalamic region in mice. Experiments of Fezf1-deficient mice indicate that Fezf1 is required for the penetration of olfactory axons through the CNS basal lamina before they innervate the OB.

抗原: Synthetic peptide of the C terminal part of Mouse Fezf1 (CGPLPSGLHQGHQ)

配方: Lyophilized product from 1 % BSA in PBS containing 0.05 % Na₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Reacts with mouse Fezf1 and not react with Fezf2.