

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

FGFR2/K-sam Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01154

第三版

描述:The K-sam gene was cloned and identified as a fibroblast growth factor receptor (FGFR) 2 gene. K-sam generates several variant transcripts by alternative splicing, including those coding for soluble forms of the receptor.

抗原:Synthetic peptide of a part of Human FGFR2/K-sam

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:React with Scirrhou carcinoma.