
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-FHIT (F130) Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01156

第三版

描述:A candidate tumor suppressor gene, FHIT, was indicated at chromosome 3p14.2 spanning the FRA3B common fragile site. Abnormalities in structure and expression of the FHIT gene have been malities in structure and expression of the FHIT gene have been detected in a considerable fraction of lung tumors of small and non-small cell types.

抗原:Synthetic peptide for C-terminal of Human FHIT

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: