

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Flotillin-1 (C) Anti-Murine Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01157

第三版

描述: It is difficult to distinguish catch and identify Raft protein, unlike Caveole which is the another micro domain of membrane, because the protein moves continuously on a cell membrane, and because the size and the structure of the protein are labile.

抗原: Synthetic peptide of common part of C terminal of rat and mouse Flotillin-1 (SQVNHNKPLRTA)

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Reacts with Rat, Mouse and Human Flotillin-1 (Confirmed by western blotting and Immunohistochemistry)