

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Gastrin Receptor (OK-524) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01159

第三版

描述:Gastrin is one of the most important gut-related hormones contributing to gastric acid secretion. Gastrin and the gastrin receptor (GR) are coexpressed in many gastrointestinal cancer cell lines, and gastrin is regarded as an autocrine growth factor for these cells. Thus, gastrin is important in the investigation of gastrointestinal carcinogenesis and anticancer therapy. The GR (identical with the cholecystokinin B receptor; CCKBR) gene was cloned by Lee et al. and has an open reading frame encoding 447 amino acids with 7 transmembrane domains. To date, many studies have revealed the localization of GR in the human stomach. It has been of public interest to investigate the expression and potential role of gastrin and GR in gastric carcinogenesis and their relation to Helicobacter pylori infection and inflammation.

抗原:Synthetic peptide of the C terminal part of Human Gastrin Receptor (CTPSIASLSRSLSYTTIS)

配方:Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: