
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-GFAT1 (663) Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01165

第三版

描述:GFAT1 (Glutamine : fructose-6-phosphate amidotransferase 1) is one of nutrientsensitive phosphoprotein which was identified as a protein phosphorylated in glucose-deprived cells.

抗原:Synthetic peptide of the C terminal part of GFAT1 (RGYDVDFPRNLAKSVTVE)

配方:Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Cross reacts with mouse and rat.