

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Gipie Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01169

第三版

描述:A new Girdin family molecule, Gipie is expressed in endothelial cells and it is induced by ER stresses. It interacts with GRP78 that is a 78kDa glucose regulatory protein. It has been reported that Gipie controls IRE1-JNK signal pathway by the interactive function and protects endothelial cells from apoptosis that is induced by ER stresses under circumstances such as atherosclerosis or vascular endothelial failure.

抗原:Synthetic peptide of Gipie (VQEPEKRPLTPSLSQ).

配方:Lyophilized product in PBS containing 1 % BSA and 0.05 % NaN₃.

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Cross-react with Human Gipie.