

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

GLUT-1 Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01173

---

第三版

**描述:** Mammalian cells transport glucose through a family of membrane proteins known as glucose transporters (GLUTs or SLC2A family). Mammalian tissues are known to express different glucose transporter isoforms (existing at least GLUTs 1-13) in different tissues. The molecular weight are 40-60 kDa, respectively.

**抗原:** Synthetic peptide of the C terminal part of human GLUT-1 (SDKTPEELFHPLGADSQV)

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Confirmed by western blotting using transfectant cells.