
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

GLUT-3 Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01175

第三版

描述:Mammalian cells transport glucose through a family of membrane proteins known as glucose transporters (GLUTs or SLC2A family). Mammalian tissues are known to express different glucose transporter isoforms (existing at least GLUTs 1-13) in different tissues. The molecular weight are 40-60 kDa, respectively.

抗原:Synthetic peptide of the C terminal part of human GLUT-3 ((C)SGKDGVMEMNSIEPAKETT)

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Confirmed by western blotting using transfectant cells.