

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

GRO/CINC-3 Anti-Rat Rabbit IgG

多克隆抗体

产品货号: PA01187

第三版

**描述:** Cytokine-induced neutrophil chemo attractant 1 (CINC-1) was originally purified from media conditioned by IL-1 stimulated rat kidney epithelioid cells (NRK-52E).

**抗原:** Recombinant Rat GRO/CINC-3.

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Rat GRO/CINC-3 specific. Non-cross react with Rat GRO/CINC-1, Rat GRO/CINC-2 and Rat GRO/CINC-2 .