



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

GRO /MGSA (C) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01188

第三版

描述:GRO /MGSA as IL-8 family has the high homology with IL-8 in amino acid sequences. It is reported that it also has the similar configuration in 3D structural analysis and the almost same level of leukocyte chemotaxis with IL-8.

抗原:Synthetic peptide in portion of C-terminus of Human GRO /MGSA.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Non-Crossreacts to Human GRO ， Human IL-8, Human IP-10, Rat GRO/CINC-1, Rat GRO/CINC-2 ， Rat GRO/CINC-2 and Rat GRO/CINC-3.