

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-GZF1 Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01192

第三版

描述:GZF1 has been identified as a gene whose expression is induced by a neurotrophic factor, GDNF. It codes a protein which has BTB/POZ domain in its N-terminal region and 10 of zinc finger motifs in its C-terminal area from the center. It has been known by performance analysis that GZF1 is a transcriptional repression factor which sequence-specifically binds to DNA and the binding sequence is present in transcription control region of HOXA10.

抗原:Synthetic peptide of the part of Human GZF1

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN3

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:React with human and mouse GZF1