
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Syndecan-4 Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01196

第三版

描述:Syndecan, a cell surface proteoglycan, is an integral membrane protein acting as a receptor for the extracellular matrix. The syndecan family of four transmembrane heparan sulfate proteoglycans binds a variety of soluble and insoluble extracellular effectors.

抗原:Maltose Binding Protein (MBP) - Syndecan-4 recombinant protein

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Confirmed by western blotting