
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

-Synuclein (S122) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01197

第三版

描述: The pre-synaptic protein, α -synuclein, has been associated with the pathogenesis of several neurodegenerative diseases. These include prevalent conditions such as Parkinson's disease (PD), dementia with Lewy bodies (DLB), and the Lewy body variant of Alzheimer's disease (LBVAD), as well as rarer conditions including multiple systems atrophy (MSA), and neurodegeneration with brain iron accumulation type-1 (NBIA-1).

抗原: Synthetic peptide for C-terminal of Human α -Synuclein (MKDEEYEQM)

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Confirmed by western blotting