

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Synuclein Antibody

多克隆抗体

产品货号: PA01199

第三版

**描述:** Synucleins are small, soluble proteins expressed primarily in neural tissue and in certain tumors. The family includes three known proteins:  $\alpha$ -synuclein,  $\beta$ -synuclein, and  $\gamma$ -synuclein. All synucleins have in common a highly conserved  $\alpha$ -helical lipid-binding motif with similarity to the class-A2 lipid-binding domains of the exchangeable apolipoproteins

**抗原:** Recombinant human  $\alpha$ -Synuclein purified from E. coli

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**如何使用:** 加1ml超纯水重溶

**稳定性:**  $-20^{\circ}\text{C}$ 保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) , 1% BSA

**应用:** WB 1 ~ 5  $\mu\text{g/ml}$ .

**特异性:** This antibody is specific for  $\alpha$ -synuclein and does not cross-react with  $\beta$ - and  $\gamma$ -synuclein.