

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Tau Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01200

第三版

**描述:** Tau proteins are microtubule-associated proteins that are thought to regulate different properties of neuronal and possibly also glial microtubules, especially their stability and orientation.

**抗原:** Synthetic peptide for N terminus portion (AEPRQEFVEMDHAG) of human Tau Protein

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Reacts with all of recombinant human Tau protein isoforms, Tau 352, 383, 381, 412, 410, and 441, which are confirmed by western blotting.