



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Tie (N1125) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01202

第三版

**描述:**This antibody can be stained in formalin fixed paraffin embedded tissues after microwave pretreatment by several Immunohistochemical techniques such as Avidin Bition Complex (ABC) method.

**抗原:**Part of synthetic peptide of Human Tie-1

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**Confirmed by Tie/BaF3 Transfectant and HEL