
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

TIMP-1 (T172) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01204

第三版

描述: Tissue inhibitor of Metalloproteinase (TIMP-1) is endogenous-specific inhibitor produced simultaneously with matrix metalloproteinase (MMP). TIMP-1 inhibits both active MMP to control degeneration of extracellular matrix (ECM) and conversion of intermediate form into active form.

抗原: The part of Synthetic peptide of TIMP-1 protein (177-188)

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Confirmed by Recombinant TIMP-1