

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Tob (Phosphorylated) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01206

---

第三版

**描述:**Tob is a member of an emerging family of genes with anti-proliferative function. Tob is rapidly phosphorylated at Ser 152, Ser 154, and Ser 164 by Erk1 and Erk2 upon growth-factor stimulation. Tob inhibits cell growth by suppressing cyclin D1 expression, which is canceled by Erk1- and Erk2-mediated Tob phosphorylation (1).

**抗原:**Synthetic peptide of human Tob protein (ASSVSS(S)P(S)PPFGHS, (S): phosphorylated)

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**Confirmed by western blotting using transfectant cells