

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

TPO Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01207

---

第三版

**描述:**Thrombopoietin (TPO) is produced primarily in the liver and secondarily in the kidneys and bone marrow. It can bind to a receptor (which is the proto-oncogene c-Mpl, a homolog of an envelope protein of the myeloproliferative leukemia virus) on primitive blood cells and megakaryocytes and stimulate those cells to grow (proliferate) and evolve into platelet-producing megakaryocytes.

**抗原:**Synthetic peptide for N-terminal of Human TPO

**配方:**Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4)，1% BSA

**应用:**WB 1 ~ 5 μg/ml.

**特异性:**