
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

TSC1/Hamartin (T1149) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01209

第三版

描述:TSC1 is a tumor-suppressor gene encoding a protein TSC1/Hamartin of about 130 kDa (ref. 1). Loss of its function causes TSC (tuberous sclerosis complex) or angiomyolipoma in human and in an animal model, TSC1 knockout mouse develops a renal tumor (ref. 2).

抗原:The synthetic peptide of a part of Human TSC1/Hamartin (GQLHIMDYNETHHEHS)

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Cross reacts with mouse.