

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

VEGF (V-N) Anti-Rat Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01217

---

第三版

**描述:** Vascular Endothelial Cell Growth Factor (VEGF) is a homodimeric protein initially purified from media conditioned by normal bovine pituitary folliculo-stellate cells and secreted by a variety of vascularized tissues.

**抗原:** Synthetic peptide for N - terminal of[ Rat VEGF

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Rat specific. Cross react with Human VEGF (165). React with rat VEGFs both of reduced form and of non-reduced form.