

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

XIAP Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01222

---

第三版

**描述:** IAP (inhibitor of apoptosis protein) is an evolutionarily-conserved inhibitory factor of apoptosis. IAP directly binds to caspase which plays the central role in operation of apoptosis and inhibits it.

**抗原:** Synthetic peptide of a part of Human XIAP (aa244-263: VSSDRNFPNSTNLPRNPSMA)

**配方:** Lyophilized product from 1 % BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** Cross reacts with mouse.