

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

GRO/CINC-2 , 2 (N) Anti-Rat Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01232

第三版

描述: Cytokine-induced neutrophil chemo attractant 1 (CINC-1) was originally purified from media conditioned by IL-1 stimulated rat kidney epithelioid cells (NRK-52E).

抗原: Synthetic peptide for N-Terminal common sequences of rat GRO/CINC-2 and rat GRO/CINC-2 conjugating bovine thyroglobulin.

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: React with Rat GRO/CINC-2 and -2 . Non-cross react with Rat GRO/CINC-1 and Rat GRO/CINC-3.