
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

HGF (H56) Anti-Rat Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01247

第三版

描述: Hepatocyte Growth factor (HGF) was discovered as a mitogen for hepatocytes. HGF was subsequently found to be identical to the scatter factor, which destroys epithelial cell adhesion and promotes cell movement.

抗原: Synthetic peptide for N-Terminal of rat HGF

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Rat HGF- specific. Non-cross react with Human HGF and Rat HGF .