



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

hnRNP B1 Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01252

第三版

描述:hnRNP proteins are a family at least 30 polypeptides, play a role in the splicing of pre-mRNA into the mature product in association with small nuclear ribonucleoprotein (snRNP) and splicing factors.

抗原:Synthetic peptide of part of N terminal of human hnRNP (KTLETVPLERKKREK)

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Confirmed by western blotting (37 kD).