

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

IL-6 (N) Anti-Rat Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01262

第三版

描述: Interleukin-6 (IL-6) is one of cytokines which was found as a B cell differentiation-inducing factor. In these days, IL-6 is generating a lot of attention not only from antibody-producing system but also from assorted physiologically activity such as induction in biosynthesis of acute stage protein in hepatocytes and acceleration of proliferation based on synergism with IL-3.

抗原: Synthetic peptide in portion of N-terminus of Rat IL-6

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: