

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

IL-8 (C) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01264

第三版

**描述:** Interleukin-8 (IL-8) is leukocyte chemotactic factor produced by various cells such as leukemia, fibroblast and endothelial cells followed by stimulation of inflammatory cytokine such as IL-1 and TNF.

**抗原:** Synthetic peptide in portion of C-terminus of Human IL-8

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:** NonCross-reacts with Human GRO Human GRO /MGSA, Rat GRO /CINC-1, 2 , 2 , 3