

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

IRF-3 Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01267

第三版

描述: Interleukin Regulatory Factor-3 (IRF-3) is the gene coding for a regulatory factor that mediates the action of IL-1, TNFalpha and LPS. It is reported that IRF-3 controls the production of interferons and chemokines in virus infected cells, by a posttranslational mechanism that involves virus induced phosphorylation, protein dimerization, cytoplasmic to nuclear translocation and stimulation of responsive genes.

抗原: Synthetic peptide for a.a.137 – a.a.150 (EDILDELLGNMVL A) of IRF-3

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Confirmed by western blotting using His tag IRF-3.