

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

IRS (Isoleucyl-tRNA synthetase)

多克隆抗体

产品货号: PA01269

---

第三版

**描述:** IRS is a component of multi-ARS complex and catalyzes the ligation of isoleucine with isoleucyl-tRNA.

**抗原:** S1074-F1264 denatured recombinant human isoleucyl tRNA synthetase protein.

**配方:** Anti-isoleucyl-tRNA synthetase rabbit serum containing 0.05% sodium azide for preservative.

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** This antibody reacts specifically with 135kDa isoleucyl tRNA synthetase protein.