

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

I B- (I037) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01270

第三版

描述:In signal transduction system, I B acts as inhibitor to forms complex with NF- B in cytoplasm and subsequently to inactivate NF- B. NF- B is transcription factor distributed in various cells.

抗原:Synthetic peptide in portion of N-terminus of Human I B- protein

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Confirmed by Human I B-