



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Kappa Light Chain Goat IgG Biotin (Monkey Absorbed)

多克隆抗体

产品货号: PA01273

第三版

描述:Anti-Human IgG (goat IgG) is often used as a secondary antibody for the detection of Human IgG. Following recent advances in the development of antibody-based therapy, such as the development of Human antibodies, this secondary antibody has been used for the purpose of investigating metabolism in the living body following administration of this antibody.

抗原:Purified Human Kappa light chain

配方:Lyophilized product from 1 % BSA in PBS containing 0.05 % Na₃N

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: