

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

KC (C) Anti-Mouse Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01274

第三版

描述: Mouse KC (Keratinocyte-derived Cytokine) is a counterpart of Human GRO protein. It is reported that KC is produced in such as activated monocyte and fibroblast.

抗原: Synthetic peptide in portion of C-terminus of Mouse KC

配方: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: