

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

KC (N) Anti-Mouse Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01275

第三版

**描述:** Mouse KC (Keratinocyto-derived Cytokine) is a counterpart of Human GRO protein. It is reported that KC is produced in such as activated monocyte and fibroblast.

**抗原:** Synthetic peptide in portion of N-terminus of Mouse KC

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μg/ml.

**特异性:**