

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

KC Anti-Mouse Rabbit IgG

多克隆抗体

产品货号: PA01277

第三版

描述:Mouse KC (Keratinocyto-derived Cytokine) is a counterpart of human GRO protein. It is reported that KC is produced in such as activated monocyte and fibroblast.

抗原:Recombinant Mouse KC

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性: