



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

KDR (C) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01279

第三版

描述: VEGF is a pleiotrophic factor that exerts a multitude of biological effects through its interaction with two receptor tyrosine kinases, Flt-1 and KDR/Flk-1. The KDR/Flk-1 is responsible for the proliferative activities of VEGF, considerable controversy and uncertainty exists about the role of the individual receptors in eliciting many of the other effects.

抗原: Synthetic peptide for C-terminal of Human KDR

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: