
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Leptin (C) Anti-Human Rabbit IgG Affinity Purify

多克隆抗体

产品货号: PA01282

第三版

描述:Leptin is a 16 KDa secreted protein produced by the ob gene. Adipose tissue produces Leptin and releases it into the bloodstream. As fat deposits grow, blood Leptin levels tends to increase.

抗原:Synthetic peptide in portion of C-terminus of Human Leptin

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Cross-reacts with Mouse and Rat