

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

LRG (128) Anti-Rat Rabbit Affinity Purify

多克隆抗体

产品货号: PA01291

第三版

描述: Among the diseases which should be distinguished from iNPH (idiopathic normal pressure hydrocephalus) of elderly people showing symptoms such as gait disorder and dementia, there are some neurodegenerative diseases such as Alzheimer's disease, FTLD (frontotemporal lobar degeneration) and disorder associated with Parkinson's disease. iNPH can be treated by ventriculo-peritoneal shunt (VP shunt), so a marker in spinal fluid which can distinguish it from neurodegenerative diseases in an early stage has been required.

抗原: Synthetic peptide of a part of Rat LRG (LVLRENQLQEASARWLQG)

配方: Lyophilized product from 1 % BSA in PBS containing 0.05% NaN₃

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Rat LRG specific.